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Guidance and manufacturer's declaration- electromagnetic emissions		
The medical suction is intended for use in the electromagnetic environment specified below. The customer or the user of the medical suction should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The medical suction uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	The medical suction is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity			
The medical suction is intended for use in the electromagnetic environment specified below. The customer or the user of the medical suction should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kv contact ±8 kv air	±6 kv contact ±8 kv air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kv for power Supply lines ±1kv for input /output lines	±2 kv for power Supply lines ±1kv for input /output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kv differential mode ±2kv common mode	±1 kv differential mode ±2kv common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) For 0.5 cycle 40% UT (60% dip in UT) For 5 cycles 70% UT (30% dip in UT) For 25 cycles <95% dip in UT) For 5 sec	<5% UT (>95% dip in UT) For 0.5 cycle 40% UT (60% dip in UT) For 5 cycles 70% UT (30% dip in UT) For 25 cycles <95% dip in UT) For 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the medical suction requires continued operation during power mains interruptions, it is recommended that the medical suction be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8			Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The medical suction is intended for use in the electromagnetic environment specified below. The customer or the user of the medical suction should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150kHz to 80 MHz</p> <p>3V/m 80 MHz to 2.5GHZ</p>	<p>[3]V</p> <p>[3]V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the medical suction, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[ \frac{3.5}{3} \right] \sqrt{120} = 12.78$ $d = \left[ \frac{3.5}{3} \right] \sqrt{120} = 12.78$ $d = \left[ \frac{7}{3} \right] \sqrt{120} = 25.56$ <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (<math>W</math>) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (<math>m</math>).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>a</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situation. Electromagnetic propagation is affected by absorption and reflection from structures, objects people.			
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic site survey should be considered. If the measured field strength in the location in which the medical suction is used exceeds the applicable RF compliance level above, The medical suction should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the medical suction.			
Over the frequency range 150kHz to 80 MHz, field strengths should be less than 10 V/m.			

Recommended separation distances between Portable and mobile RF communications equipment and the medical suction			
The medical suction is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer of the user of the medical suction can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical suction as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{V1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E1} \right] \sqrt{P}$
0.01	0.1166	0.1166	0.2333
0.1	0.3689	0.3689	0.7379
1	1.167	1.167	2.333
10	3.689	3.689	7.3786
100	11.667	11.667	23.333
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter. Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacture.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			